

This headrest is designed to support the head and neck of the person being lifted. It may be hard to keep the head upright whilst reassuring the fallen person. Therefore, we recommend this headrest to provide complete support. It is extremely useful in cases where it is necessary or practical for the helper to have both hands free.

Standards

The headrest is in compliance with relevant requirements in the following standards:

- DIN EN ISO 10993-10:2014-10 Biological Evaluation Of Medical Devices - Part 10: Tests For Irritation And Skin Sensitization
- DIN EN ISO 10993-5:2009-10 Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity











How to use the headrest





Liftup reserves the right to make technical modifications.

JD HEALTHCARE Group

To place a purchase order, contact us today: sales@jdhealthcare.com.au | 1300 791 404



